

Low speed centrifuge CD-0506

The new centrifuge CD-0506, supplied with a built in fixed angle rotor which can hold up to 6 tubes each 15ml/10ml/7ml/5ml/ 1,5/2ml. It is widely used for blood or urine separation in medical routine practice in hospitals, clinics or research institutes.



CD-0506

All features at a glance

- Speed range from 300-5000rpm
- Rotor speed can be set and displayed by RPM or G-force
- Defined program by knob will be stored and activated when power on
- Two programs P1/P2, easy to start the procedure by one key
- Two deceleration speeds setting to meet different requirements
- Easy-to-read processing display and sound alert
- Automatic lid-lock release after running
- Maintenance-free brushless DC motor
- Rotor and Adapter are included
- 2 years warranty



Adapter
CD05-A2 and CD05-A5

Technical Data	CD-0506
Speed	300-5000 rpm
Increment	10 rpm
max. RCF	2.600 g, increment 10 rpm
Timer	0,5-99 minutes, continuous
max. Capacity	6 x 15 ml
Safety Devices	Door interlock, overspeed detector, automatic internal diagnosis
Dimension	24 x 30 x 18 cm (D x W x H)
Weight	5,2 kg
Rotor	54° fix angel
Noise	< 54 dB (A)
Voltage	100 - 240 V
Frequency	50/60 Hz



Issue summer 2021, subject to alterations

