PHOENIX Instrument
Lab Solutions

ASN / ABN - Series

The ASN / ABN-Series is the all-round analytical balance for the laboratory. It guarantees brilliant stability, response and reliability characteristics. The innovative and compact construction of the weighing cell in relation with the digital process technology enables fastest weighing results and highest stability.



ASN-224

ASN / ABN-Series

- huge, user-friendly weighing chamber
- easy handling of all procedures
- RS-232 interface
- direct data transfer into Windows application by Windows direct function
- various application programs (e. g. formulation, counting, %)
- GLP/GMP-conform protocol
- user-defined printing intervals
- aluminum die-cast housing
- protection cover as standard

3 Years Warranty

Model	Accuracy	Capacity	Pan Size	Calibration	Meas. Time
ASN-324	0.1 mg	320 g	90 mm	internal	2 - 3 sec
ASN-224	0.1 mg	220 g	90 mm	internal	3 sec
ASN-124	0.1 mg	120 g	90 mm	internal	3 sec
ABN-224	0.1 mg	220 g	90 mm	extern	3 sec
ABN-124	0.1 mg	120 g	90 mm	extern	3 sec

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de