PHOENIX Instrument
Lab Solutions

Moisture Analyzer BM-60



BM-60

BM-60

- the large, backlit graphic display for easily read out even under poor lighting conditions
- the huge sample pan allows placing even larger samples on the instrument
- infrared quartz heater provides effective and homogen drying without interference
- auto-zero function ensures high accuracy measurements
- graphic display of the drying curves
- 5 different drying modes can be combined with 2 stop-functions



Technical Data	BM-60
Capacity	60 g
Accuracy	0.001 g / 0.01 %
Reproducibility	0.002 g / 0.05 %
Pan size	95 mm
Meas. units	Gram, % Moisture, % Solids
Stabilisation time	2-3 sec
Internal memory	10 Programs, 100 Test results
Calibration	automatic external
Heating source	400 W-Halogen heater, 50 - 160 °C
Drying programs	Standard, Drying in steps, gentle and fast huge
Display	large backlit LCD display

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de