PHOENIX Instrument
Lab Solutions

Moisture Analyzer BM-65

The newly developed moisture analyzer BM-65 is characterized by numerous features, helping the user to make working and testing with the instrument easy and convenient.

BM-65

- large, backlit display for easy readouts
- various parameters can be displayed simultaneously
- automatic test function
- powerful 400 W halogen heater
- internal memory
- automatic test run after the lid is closed
- pan lifter makes placing and changing of probes easier
- bi-lingual menue
- RS232 interface



BM-65



Technical Data	BM-65
Capacity	65 g
Accuracy	0.001 g / 0.01 %
Reproducibility	0.002 g / 0.05 %
Pan size	90 mm
Meas. units	Gram, % moisture, % solids, % moisture on dry bases (ATO %S)
Stabilisation time	4 sec
Internal memory	Up to 99 test results or 49 programable user setting
Calibration	external
Heating source	400 W-Halogen heater, 50 - 160 °C
Drying programs	standard one step up to 3 temp. stepping ramps
Display	large backlit LCD display

3 Years Warranty

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de