PHOENIX Instrument
Lab Solutions

## BTB-Series

The BTB-Series is equipped with the revolutionary UNIBLOC-weighing system which guarantees brilliant stability, response and reliability characteristics.
The innovative and compact construction of the weighing cell in relation with the digital process technology enables fastest weighing results and highest stability. Calibration per push on the CAL-button with an external calibration weight.

### BTB-Series

- automatic adjustment to environmental conditions
- easy handling of all procedures
- RS-232 interface
- weighing in different units
- weighing below the balance
- direct data transfer into Windows application by Windows direct function
- various application programs (e. g. formulation, counting, %, animal weighing ...)
- user-defined printing intervals
- aluminum die-cast housing
- backlit display
- protection cover as standard






5 Years Warranty

BTB-6202

| Model | Accuracy | Capacity | Pan Size |
|---|---|---|---|
| BTB-623 | 0.001 g | 620 g | 108 x 105 mm |
| BTB-423 | 0.001 g | 420 g | 108 x 105 mm |
| BTB-6202 | 0.01 g | 6,200 g | 170 x 180 mm |
| BTB-4202 | 0.01 g | 4,200 g | 170 x 180 mm |
| BTB-2202 | 0.01 g | 2,200 g | 170 x 180 mm |

External Calibration
3 seconds measuring time

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de