PHOENIX Instrument
Lab Solutions

BTC-Series

The brandnew BTC-Series sets new standards with the revolutionary weighing system UNIBLOC and the state-of-the-art design.

5 Years Warranty



BTC-423V



BTC-Series

- built-in calibration weight
- separate stability indication
- RS-232 interface
- huge weighing chamber with glass wind protection as standard
- fast and convenient handling with the new menu operating key
- direct data transfer into Windows application by Windows direct function
- various application programs (e. g. formulation, counting, %)
- new backlit display for more comfortable readings

Internal Calibration
2 seconds measuring time

Model	Accuracy	Capacity	Pan Size
BTC-423	0.001 g	420 g	110 mm
BTC-323	0.001 g	320 g	110 mm
BTC-223	0.001 g	220 g	110 mm

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de