PHOENIX Instrument
Lab Solutions

BTD-Series

The BTD-Series combines state-of the art design, easy handling, extensive equipment and brilliant weighing results delivered by using the revolutionary UNIBLOC-weighing system, to a unique standard balance.

BTD - Series

- calibration per push on the CAL-button with an external calibration weight
- separate stability indication
- RS 232 interface
- huge weighing chamber with glass wind protection as standard (mg-balances)
- fast and convenient handling with the new menu operating key
- direct data transfer into Windows application by Windows direct function
- various application programs (e. g. formulation, counting, %)
- new backlit display for more comfortable readings

5 Years Warranty





BTD-3202

Model	Accuracy	Capacity	Pan Size
BTD-423	0.001 g	420 g	110 mm
BTD-323	0.001 g	320 g	110 mm
BTD-223	0.001 g	220 g	110 mm
BTD-4202	0.01 g	4,200 g	167 x 181 mm
BTD-3202	0.01 g	3,200 g	167 x 181 mm
BTD-2202	0.01 g	2,200 g	167 x 181 mm

External Calibration
2 seconds measuring time

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de