PHOENIX Instrument
Lab Solutions

BTE-Series

State-of-art-design, smart handling, extensive equipment, top prices.

BTE-Series

- calibration per push on the CAL-button with an external calibration weight
- separate stability indication
- fast and convenient handling with the new menu operating key
- direct data transfer into Windows application by Windows direct function
- RS 232 interface
- various application programs (e. g. formulation, counting, %)
- new backlit display for more comfortable readings



BTE-622

External Calibration
2 seconds measuring time

Model	Accuracy	Capacity	Pan Size
BTE-622	0.01 g	620 g	110 mm
BTE-6201	0.1 g	6,200 g	160 mm
BTE-4201	0.1 g	4,200 g	160 mm

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de