PHOENIX Instrument
Lab Solutions

BTF-Series

The brandnew BTF-series offers an excellent economical alternative in the range of toploading balances with an accuracy of 0.001 g to 0.01 g.

BTF-Series

- internal calibration weight (BTF-1003i, BTF6202i, BTF-4202i)
- huge, backlit display
- glass ring protection as standard at mg-balances
- convenient and easy operation
- various application programs
- RS-232-interface

3 Years Warranty



BTF-1003

Model	Accuracy	Capacity	Pan Size	Calibration
BTF-1003	0.001 g	1,000 g	110 mm	extern
BTF-1003i	0.001 g	1,000 g	110 mm	inten
BTF-6202	0.01 g	6,200 g	160 mm	extern
BTF-6202i	0.01 g	6,200 g	160 mm	inten
BTF-4202	0.01 g	4,200 g	160 mm	extern
BTF-4202i	0.01 g	4,200 g	160 mm	inten

3 seconds measuring time

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de