PHOENIX Instrument
Lab Solutions

BTG-224

The new neue analytical balance BTG-224 is ideal for all lab purpose. It combines easy handling, modern design and precise results.



BTG-224

BTG-224

- internal calibration
- reliable weighing system
- GLP-conform protocol (when connected to PC or printer)
- different user programs
- RS-232 Interface
- protection cover

3 Years Warranty

Technical Data	BTG-224
Capacity	220 g
Accuracy	0.1 mg
Reproducibility	0.1 mg
Linearity	0.2 mg
Meas. time	3 sec
Rechargeable Batteries	yes
Pan size	90 mm
Dimension	yes

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de