PHOENIX Instrument
Lab Solutions

## BTG-Series

The brandnew BTG-series consists of multi-functional balances with a number of application programs.They are easy to handle and well-equipped. The internal rechargeable battery ensures a maximum of portability.

### BTG-Series

- large, bright display
- intuitive and easy operation, large weighing pan
- Model BTG-303 with protection housing
- Calibration weight included
- various application programs
- RS-232 interface

3 Years Warranty



BTG-2002







BTG-303

External Calibration

| Model | Accuracy | Capacity | Pan Size | Meas. Time |
|---|---|---|---|---|
| BTG-303 | 0.001 g | 300 g | 118 mm | 2-3 sec |
| BTG-2002 | 0.01 g | 2,000 g | 188 x 148 mm | 2 sec |
| BTG-6001 | 0.1 g | 6,000 g | 188 x 148 mm | 2 sec |

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de