PHOENIX Instrument
Lab Solutions

Entry-Level Model CB-2002 N and CB-6001 N

New entry-level balance for fast and reliable results at an affordable price.

CB-Serie

- backlit LCD - display
- rechargeable battery included
- various programs
- external calibration
- stainless steel pan

2 Years Warranty



CB-2002 N

Technical Data	CB-2002 N	CB-6001 N
Capacity	2,000 g	6,000 g
Accuracy	0.01 g	0.1 g
Reproductibility	0.02 g	0.2 g
Linearity	+/- 0.03 g	+/- 0.2 g
Meas. time	3 sec	1.5 sec
Rechargeable batteries	yes	yes
Pan size	160 x 120 mm	192 x 143 mm

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de