PHOENIX Instrument
Lab Solutions

Mini-Centrifuge CD-1008

The CD-1008 is a small and powerful mini centrifuge which fits on every laboratory bench. Compact, small and very quiet running (<45 dB) is it the ideal instrument for small analysis. For safety reasons it is running only with a closed lit.

CD-1008 G (green), **CD-1008 B** (blue)

- Two lid colors, green and blue
- Max. speed of 7,000 rpm
- Comes with 2 ml x 8 and PCR8 x 2 rotors
- Dual door interlock design
- DC motor drive
- Quiet running, noise level ≤45 dB
- The rotor is seated on a tapered shaft which facilitates easy removal
- Ideal for quick spin downs and micro filtration
- The centrifuge starts and stops by closing and opening lid or using the on/off switch



2 Years Warranty



Accessories

CD08-8-2	Rotor, 8 x 1.5/2 ml
CD08-2-PCR8	Rotor, 2 x 8 x 0.2 ml PCR
CD08-A02	Adapter, 2 ml to 0.2 ml, 8 pcs
CD08-A05	Adapter, 2 ml to 0.5 ml, 8 pcs

Technical Data	CD-1008
Speed	7,000 rpm
Max. RCF	2,680 g
Time	continous
Noise	<45 dB (A)
Dimension	15 x 15 x 12 cm
Weight	0.5 kg
Max. capacity	8 x 2 ml
Shipping includes	CD 08-2, CD 02-PCR 8
Voltage • Frequency	100 - 240 V • 50/60 Hz
Power	20 W

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de