PHOENIX Instrument
Lab Solutions

Digital Sea-Saw Shaker RS-RR 10 and 3D-Shaker RS-RD 10

RS-RR 10 • RS-RD 10

- gentle sample mixing
- sturdy construction
- backlit display
- simultaneous display of speed and time
- fix stand with anti-slipping feet
- attachment included
- comprehensive range of accessories
- RSS-E 220 included





2 years warranty

Technical Data	RS-RR 10	RS-RD 10
Shaking movement	Sea-saw	3D orbital
Shaking angel	±9 °	
Max. speed	10 - 70 rpm	
Dispaly	LCD	
Timer-function	yes, 1 - 1.199 min	
Operation	Timer / Continous	
Max. shaking weigt	10 kg	5 kg
Dimension	420 x 370 x 106 mm	420 x 370 x 200 mm
Dimension Platform	320 x 320 mm	
Ambient temperature	5 - 40 °C	
Protection class	IP 21	
Voltage • Frequency	100 - 240 VAC • 50/60 Hz	
Weight	13.9 kg	15.1 kg

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de