PHOENIX Instrument
Lab Solutions

Rotator RS-RR 20 and Disk Mixer RS-RD 20


2 years warranty



RS-RR 20

RS-RR 20 • RS RD 20

- instrument will be supplied incl. tube holder (RS-RR20) or disk mixer (RS-RD 20)
- dual-Use - rotator RS-RR 20 can also be used as disk mixer
- fix stand with anti-slipping feet
- backlit display to adjust speed and time
- comprehensive range of accessories



RS-RD 20

Technical Data	RS-RR 20	RS-RD 20
Movement	Rotation	Rotation Disk
Rotation angle	-	0 - 90 °
Speed	10 - 70 rpm	
Timer-function	yes	
Operation	Timer / Continous	
Dispaly Speed / Time	LCD / LCD	
Adjustment Speed / Time	digital / digital	
Dimension	510 x 190 x 260 mm	250 x 190 x 260 mm
Protection class	IP 21	
Power	40 W	
Voltage • Frequency	100 - 240 VAC • 50/60 Hz	
Weight	8.9 kg	5.8 kg

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de