PHOENIX Instrument
Lab Solutions

Analog Rotator and Disk Mixer RS-RR 5 • RS-RD 5

RS-RR 5 • RS-RD 5

- accessories for 1.5/2 ml, 15 ml and 50 ml tubes are included
- variable Speed
- compact construction
- adjustable feet
- easy to clean

2 years warranty



RS-RR 5



RS-RD 5

Maximum assembling RS-RR 5 and RS-RD 5

Fixing clips	1.5 / 2 ml	15 ml	50 ml
RS-RR 5	34	20	12
RS-RD 5	48	8	8

Technical Data	RS-RR 5	RS-RD 5
Movement	Rotation	Rotation Disk
Rotation angle	-	fix angle
Speed	5 - 80 rpm	
Timer-function	-	
Operation	continous	
Dispaly Speed / Time	Scala / Scala	
Adjustment Speed / Time	analog / analog	
Dimension	530 x 150 x 190 mm	2,200 x 300 x 310 mm
Protection class	IP 21	
Power	20 W	
Voltage • Frequency	100 - 240 VAC • 50/60 Hz	
Weight	2.9 kg	3.1 kg

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de