PHOENIX Instrument
Lab Solutions

Tube Roller RS-TR 5 / RS-TR 10



RS-TR 5

RS-TR 5 / RS-TR 10

- gentle mixing
- easy to use
- easy cleaning
- timer function (RS-TR 10 only)
- LCD display (RS-TR 10 only)



2 years warranty



RS-TR 10

Technical Data	RS-TR 5	RS-TR 10
Movement	Rocking amd rolling	
No. of rolls / Length	6 / 280 mm	
Amplitude	24 mm	
Speed	20 - 70 rpm	10 - 70 rpm
Motor	DC, brushless	
Timer-function	-	yes
Operation	continous	Timer / continous
Display Speed / Time	Scale	LCD / LCD
Adjustment Speed / Time	analog / -	digital / digital
Dimension	260 x 490 x 110 mm	
Protection class	IP 21	
Power	25 W	30 W
Voltage • Frequency	100 - 240 VAC • 50/60 Hz	
Weight	5.5 kg	5.6 kg

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de