PHOENIX Instrument
Lab Solutions

Vortexer RS-VF 10 • RS VA 10 and Micro Plate Mixer RS-MM 10



RS-MM 10

- compact Mixer suitable for 96er Micro titer plates
- adjustable speed
- various top frames available for RS-MM 10

2 years warranty



RS-VF 10 • RS-VA 10

- touch- or continuous
- variable speed 0 - 3,000 rpm (RS-VA 10)
- anti-slipping feet
- various top frames available for RS-VA 10
- sturdy cast housing

Technical Data	RS-VF 10	RS-VA 10	RS-MM 10
Shaking movement	orbital		
Amplitude	4 mm		4.5 mm
Speed	3,000 rpm	100 - 3,000 rpm	50 - 1,500 rpm
Speed adjustment	-	knob	knob
Max. load	-		0.5 kg
Dimension	127 x 160 x 130 mm		150 x 80 x 260 mm
Protection class acc. EN60529	IP 21		
Dimension	155 x 215 x 100 mm		180 x 320 x 108 mm
Weight	3.5 kg		
Voltage • Frequency	100 - 240 VAC • 50/60 Hz		

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de

PHOENIX Instrument
Lab solutions

Vortexer and Microplate Mixer - Accessories


RSV-E 10	Universal top plate for tubes and vessels up to 30 mm Ø for Vortexer
RSV-E 12	Holding bar for RSV-E 20 – 40
RSV-E 15	100 mm plate with rubber mat for RSV-E 20-40, RSV-E 44
RSV-E 18	DNA -attachment 20 Tubes 10,5 mm Ø, Material PP
RSV-E 20	Adapter for 48 Tubes 6 mm Ø, to be used with RSV-E 15
RSV-E 21	Adapter for 56 Tubes 5 mm Ø, 7 x PCR-Strips, to be used with RSV-E 15
RSV-E 25	Adapter for 15 Tubes 10 mm Ø, to be used with RSV-E 15
RSV-E 26	Adapter for 12 Tubes Ø 9,8 mm (1.5/2 ml), 12 Tubes Ø 6.6 mm (0.5 ml), 12 Tubes Ø 5 mm (0.2 ml) to be used with RSV-E 15
RSV-E 30	Adapter for 12 Tubes 12 mm Ø, to be used with RSV-E 15
RSV-E 35	Adapter for 8 Tubes 16 mm Ø, to be used with RSV-E 15
RSV-E 40	Adapter for 8 Tubes 20 mm Ø, to be used with RSV-E 15
RSV-E 41	Adapter for 4 Tubes 28 mm Ø (50 ml), to be used with RSV-E 15
RSV-E 42	Tube-Adapter horizontal, 4 x 0.5 - 15 ml
RSV-E 43	Attachment for Erlenmeyer flasks 100 - 250 ml
RSV-E 44	Universal top plate for tubes and vessels up to 99 mm Ø for Vortexer
RSV-E 46	Rubber foot
RSV-E 50	Attachment for micro plate mixer to be used with RSV-E20-40, RSV-E44
RSV-E 55	Microtiterplate-attachment for one plate for RS-MM 10
RSV-E 57	Microtiterplate-attachment for two plates for RS-MM 10
RSV-E 60	Vortexer universal attachment for tubes and vessels up to 99 mm Ø incl. RSV-E 15
RSV-E 65	Vortexer attachment with Microtiterplate RSV-E55



RSV-E 15



RSV-E 20-40



RSV-E 42



RSV-E 43



RSV-E 44



RSV-E 65

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de