PHOENIX Instrument
Lab Solutions

# Water Baths WB-5 / WB-12 / WB-22

from RT+5 to 99 °C

The waterbaths are generally used in the laboratory to maintain the temperature of the samples constant. The water-baths WB series, thanks to the maximum operating temperature of 99 °C, meet the different needs of operators and therefore allow them to be used in numerous applications.

The model WB 22 pump ensures a faster and more uniform temperature distribution due to the recirculation pump which is equipped.

All water baths are equipped with a water drain.




WB-12 with rackTIN-1001822
(max. capacity 4 pcs.)




| Technical Data | WB-5 | WB-12 | WB-22 |
|---|---|---|---|
| Volume | 5 Liters | 12 Liters | 22 Liters |
| Max. Temperature/Resolution | RT+5 to 99 °C / 0.1 °C | | |
| Temp. homogeneity at 37 °C | ±0.5 °C | | |
| Temperature variation at 37 °C | ±0.1 °C | | |
| Overheating protection | yes | | |
| Timer | 99:59 / continuous operation | | |
| Bottom plate Dimension (WxD) | 270 x 120 mm | 290 x 230 mm | 470 x 290 mm |
| Max. number of Modules | 2 | 4 | 8 |
| Min. height with closed lid | 150 mm | | |
| Voltage / Power | 230 V / 900 W | | 230 V / 1,100 W |
| Ext. Dimension (WxHxD) in mm | 480 x 215 x 335 | 480 x 380 x 310 | 680 x 390 x 365 |
| Weight | 8 kg | 12 kg | 18 kg |
| Safety class | 2 | | |




Water drain




Maximum modules in place

info@phoenix-instrument.de • Tel. +49 (5131) 90818 -30
Heinkelstr. 4 • D-30827 Garbsen • www.phoenix-instrument.de